



Superior Implant Technology

# CAD CAM PROSTHETIC SOLUTIONS

AB's CAD CAM technology for customized abutments offers products and services from CAD to CAM, from implant to final crown, from temporary restoration up to highly esthetic prosthetics.


AB's custom products are manufactured from high-quality, high-performance materials.


































Designed for precise fit with original AB implants, each product meets most stringent demands in implant dentistry for fit, stability and successful long-term results.





Superior Implant Technology



	SCAN ABUTMENTS	TITANIUM BASE FOR CROWNS	TITANIUM BASE FOR BRIDGES	PRINTED MODEL ANALOG
STANDARD PLATFORM (3.75mm diameter)	   <p>P3-3.75,SC Scan Abutment</p> <p>P3L-3.75,SC Narrow Long Scan Abutment</p>	    <p>P3-3.75,TI Slim H-0.6mm</p> <p>P3-3.75,2,TI Slim H-2mm</p> <p>P3-3.75,3,TI Slim H-3mm</p> <p>P3H-3.75,TIT Wide</p>	    <p>P3-3.75,TIT Wide Rotational</p> <p>P3N-3.75,0.6,TIT Slim Rotational H0.6mm</p> <p>P3N-3.75,2,TIT Slim Rotational H2mm</p> <p>P3N-3.75,3,TIT Slim Rotational H3mm</p>	 <p>D1-3.75,MA Digital Model Analog</p>
NARROW PLATFORM (3mm diameter)	   <p>P3-3,SC Scan Abutment</p> <p>P3L-3,SC Narrow Long Scan Abutment</p>	    <p>P3-3,TI Slim H-0.8mm</p> <p>P3-3,2,TI Slim H-2mm</p> <p>P3-3,3,TI Slim H-3mm</p> <p>P3H-3,TIT Wide</p>	    <p>P3-3,TIT Wide Rotational</p> <p>P3N-3,0.6,TIT Slim Rotational H0.6mm</p> <p>P3N-3,2,TIT Slim Rotational H2mm</p> <p>P3N-3,3,TIT Slim Rotational H3mm</p>	 <p>D1-3,MA Digital Model Analog</p>
CONICAL PLATFORM	  <p>P3LC-4,SC Scan Abutment</p>	   <p>P3NC-4,1.5 Slim H-1.5mm</p> <p>P3NC-4,2.5 Slim H-2.5mm</p> <p>P3NC-4,3.5 Slim H-3.5mm</p>	   <p>P3NC-4,1.5,R Slim H-1.5mm</p> <p>P3NC-4,2.5,R Slim H-2.5mm</p> <p>P3NC-4,3.5,R Slim H-3.5mm</p>	 <p>D1C-4,MA Digital Model Analog</p>




The CAD CAM libraries were validated in the following system: 3Shape, ZirkonZhan, Exocad, DentalWings. The most updated CADCAM products libraries can be downloaded from AB Dental website.





Superior Implant Technology



MULTI UNITS	SCAN BODEIS	ADHESIVE SLEEVES FOR CROWNS	ADHESIVE SLEEVES FOR BRIDGES	PRINTED MODEL ANALOG
<p>P14/P16</p> 	<p>P14,SC Scan Body</p> 	<p>P14-bTs Titanium Straight Sleeve</p> 	<p>P14-bTs-C Titanium Conical Sleeve</p>  <p><b>NEW</b></p>	<p>D1-P14,MA Digital Model Analog</p>  <p><b>NEW</b></p>
<p>P64</p> 	<p>P64, SC Scan Body</p> 	<p>P64-bTs Titanium Straight Sleeve</p> 	<p>P64-bTs-C Titanium Conical Sleeve</p> 	<p>D1-P64,MA Digital Model Analog</p> 
<p>P14C/P16C</p> 	<p>P14C-AR,SC Scan Body <b>NEW</b></p> <p>P14C,SC Scan Body <b>UPDATED</b></p> 	<p>P14C-bTs-AR Titanium Anti-Rotational Sleeve</p> 	<p>P14C-bTs-C Titanium Conical Sleeve</p> 	<p>D1-P14C,MA Digital Model Analog</p> 



The CAD CAM libraries were validated in the following system: 3Shape, ZirkonZhan, Exocad, DentalWings. The most updated CAD/CAM products libraries can be downloaded from AB Dental website.

